

POINT-OF-CARE

aptitude
Metrix®

Provided by
SEKISUI
DIAGNOSTICS

Metrix® COVID-19 Test

A MOLECULAR LAB ANYWHERE AND EVERYWHERE

Features and Benefits

- **Rapid Molecular Test** - Provides lab quality results in 30 minutes
- **Reduce repeat testing** - No need for serial or confirmation testing*
- **Convenient** - Indicated for symptomatic and asymptomatic testing and requires no maintenance or calibration, ready for testing any time
- **Minimal Hands-On Time** - Can be performed in a CLIA-waived setting in a few simple steps and stored at room temperature
- **Flexible sampling** - The only test with dual-sample types (nasal swab or saliva)
- **Exceptional Support** - Experienced support teams and online training modules to help streamline your implementation

Performance Characteristics

No serial testing required. No confirmatory testing needed*. A result you can trust the first time.

Sample Type	Positive Percent Agreement** (PPA)	Negative Percent Agreement** (NPA)	Limit of Detection***	Time To Results
Nasal Swab	97%	99%	2000 GE/Swab (667 GE/mL VTM equivalent)	30 Minutes
Saliva	90%	99%	20 GE/uL saliva sample	30 Minutes

*Nasal samples only. **Rounded numbers compared to RT-PCR. See IFU. ***Refer to IFU for additional LOD information.



PPA

NPA

NASAL

97%

99%

SALIVA

90%

99%



This product has not been FDA cleared or approved, but has been authorized for emergency use by FDA under an EUA; This product has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens. The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b) (1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb 3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

Metrix® COVID-19 Test

Experience the Simplicity of Metrix



1. SAMPLE

Collect nasal swab or saliva sample in the collector.



2. SNAP

Snap the cap onto the collector and then onto the sensor.



3. SLIDE

Slide the sensor onto the reader.
Read result within 30 minutes.

Interpretation of Results

Multicolor test status indicators with direct printed labels, removes human error and judgment.



Ordering Information

Test Description	Configuration	Catalog Number	CPT Code
Metrix Reader (Gen 2)	1 Reader	MTRX-RDR-Gen2	
Metrix COVID-19 Test	25 Tests	MTRX-C19-25PK	87635QW
Metrix COVID-19 Molecular Bundle	1 Reader + 50 Tests	MTRX-C19-B1	
Microbix COVID-19 Controls (for use with Metrix COVID-19 Test)	1 SARS-CoV-2 Positive Swab, 1 SARS-CoV-2 Negative Swab	RED-S-19/99-P-01	

To learn more about the Metrix Test,
scan the code or go to metrixmolecular.com



Provided by
SEKISUI Diagnostics, LLC
One Wall Street
Burlington, MA 01803
Phone: 800 332 1042
Email: questions@sekisui-dx.com
sekisuidiagnostics.com

Designed and manufactured by
Aptitude Medical Systems Inc.
125 Cremona Drive, Suite 100
Goleta, CA 93117



80-8958-00-03 04/25