



STAT-Site® WB

Glucose and β -ketone measurement

CLIA waived

Handheld & Portable analyzer

Results in just five-ten seconds

Stores 400 results

Seamless data download with EKF Link

Scan here to view



product webpage

EKF
Diagnostics

Simple and intuitive

- Auto start when the strip is inserted
- Auto reaction time countdown
- Auto switch-off after three minutes without action
- Touch button to eject test strip
- Two 1.5V AAA alkaline batteries
- 400 results stored on the system

Accurate with excellent linear range

- **Measurement range:**
 β -ketone: 0.1 to 8.0 mmol/L
 Blood glucose: 10 to 600 mg/dL (0.5 to 33.3 mmol/L)
- **Hematocrit range:**
 β -ketone: 10 to 70%
 Blood glucose: 20 to 70%
- Calibrated using the gold standard Liquicolor beta-Hydroxybutyrate method

Fast and efficient sample handling

- 1 μ L sample size capillary or venous whole blood
- β -ketone result in ten seconds
- Glucose result in five seconds

Ergonomic design

- Compact: 95.5 (l) x 57 (w) x 15 (h) mm, 58 g
- Simple, large LCD screen interface

Catalog no.	Name
7000-001	• STAT-Site® WB Dual Analyte Measurement
7010-050	• Blood Glucose Strip box 50pcs
7020-025	• β -Ketone Strip box 25pcs
7020-050	• β -Ketone box 50pcs
7014-204	• Blood Glucose Control Solution 1 x W2 (norm), 1 x B3 (high)
7024-204	• β -Ketone Control solution 1 x L1, 1 x L2



1 Insert test strip into the meter.



2 Add blood sample.



3 Result within 10 seconds for β -ketone and 5 seconds for glucose.



Tel: 001 (830) 249 0772
 Email: sales2@ekfdiagnostics.com
 ekfdiagnostics.com

Distributed by